

## INCUBATOR



### Key specification:-

Capacity	1 cu ft	1.6 cu ft	1.6 cu ft	3.3cu ft
Volume	30Ltrs	45Ltrs	95Ltrs	125Ltrs
Temp. range	5°C to 90°C			
Temp. controller	PID controlled			
Display	LED digital			
Outer body	MS powder coated			
Inner body	SS 304			
Safety feature	Over temp. thermostat			

Power supply	220v , AC ,single phase
Shelves	2 to 5 as per requirement
Optional	<ul style="list-style-type: none"><li>- Castor wheels</li><li>- Digital timer</li><li>- Fluorescent lamp</li></ul>

[www.oracleequipments.com](http://www.oracleequipments.com)