

BOD Incubator



Key specifications:-

Volume	68L	113L	173L	238L
Capacity	2.4 cu ft	4 cu ft	6.1 cu ft	10 cu ft
External chamber	SS 304/316			
Internal chamber	SS 304			
Temp. range	2°C to 60°C			
Display	LED digital			
Temp. control	PID controlled			
Sensor	PT100 sensor			

Shelves / racks	2 to 5 as per requirement (removable)
Power	220v , AC , single phase
Equip. movement	Castor wheels 4no.
Temp. accuracy	±0.5°c

www.oracleequipments.com