

## **HOT AIR OVEN**



### **Key Specifications:-**

Inner Chamber Size	18x18x18, 24x24x24, as per requirement
Display	Led display
Accuracy	±1°
Least count/Resolution	0.1°c
Power	220v , AC, Single phase
Temperature Range	50-300°c as per requirement
Air circulating Fan	1no.

### **Highlights:-**

No.of Trays	2 no.
Inner Body	GI
Outer Body	MS Powder Coated
Finish	Blue Powder Coated

## **Optional Features:-**

- 1) Extra Trays
- 2) Pair of gloves

[www.oracleequipments.com](http://www.oracleequipments.com)