

MELT FLOW INDEX TESTER



Key specifications: -

Display	LED Digital
Power	220v ,AC single phase
Temp. Control	PID controller
Temp. Range	Ambient to 300 ^o c, ambient to 400 ^o c
Accuracy	±1 ^o c
Least Count	0.1 ^o c
Sample Cutting	Automatic
Timer	yes
Material of Heating Barrel	Brass

Highlights: -

No. of Dead Weights	3no.
Body Material	Mild steel

Outer Dimensions	
Finish	Blue Powder Coated

Optional Features: -

- 1- Other specific dead weights.
- 2- HMI model of MFI.
- 3- Weighing balance of least count 0.01g.

www.oracleequipments.com