

## **MOR (Modulus Of Rupture)**



### **Key specifications:-**

Display	LED digital
Maximum capacity	500 to 700 kg
Tiles sizes	200x200mm, 300x300mm, 600x600mm, 800x800mm, 1200x1200mm & various other tile sizes
Power	220v, AC, single phase

Least count	1 Newton
Accuracy	0.2% of full scale
Load unit	Newton
Auto-stop	Yes , after tile failure (break)
Body	MS powder coated

[www.oracleequipments.com](http://www.oracleequipments.com)