

MUFFLE FURNACE



Key Specifications:-

Display	LED Digital
Inner chamber sizes	4x4x9, 6x6x12
Power	220v , AC single phase
Temp. Control	PID controlled
Temp. Range	60-1000 ^o c
Accuracy	±2 ^o c
Least count resolution	1 ^o c

Highlights: -

Insulation	Mineral glass wool
Body Material	Mild steel Powder coated
Outer Dimensions	
Finish	Blue power coated
Exhaust Port	Yes

Optional Features: -

- 1- Tongs
- 2- Gloves
- 3- Two crucibles

www.oracleequipments.com