

## **WATER BATH**

### **Key specifications:-**

Display	LED Digital
Temp. Range	30°C to 70°C
Temp. Accuracy	±1°C
Temp. Least Count	0.1°C
Power	220v, AC, Single phase
Inner Chamber size	800x800mm, 600x600mm, 800x1200mm, 800x1600mm, 1200x1600mm and various other sizes
Inner Chamber Material	SS 304
Outer Body Material	Mild Steel
Water Circulation	Yes
Castor Wheels	Yes , no.4
No. of samples tested in one go	No.5