

## **VERTICAL AUTOCLAVE**

### **(Digital)**



### **Key Specifications:-**

Internal tank capacity	22L	40L	50L	95L
Load	2kw	3kw	3kw	4kw
Chamber dimensions	10"x18"	12"x22"	14"x22"	18"x24"
No. of drums & sizes		2(9"x11")	2(9"x11")	2(12"x15")
Working pressure	15psi - 17psi( can be upgraded to 30 psi )			
Sterilizing temperature	121°c to 134°c (150°c optional)			
Sterilization time	25-30 mins			

Display	LED display
Digital timer	0-99 mins
Construction	Double/Triple walled
Outer Body	SS 304
Insulation	Mineral glass wool
Door Locking type	Radial type
Safety features	Safety valve , Radial locking, low water detector

[www.oracleequipments.com](http://www.oracleequipments.com)