

BIOSAFETY CABINET



Key specifications:-

Working area	2x2x2 ft.	3x2x2 ft.	4x2x2 ft.	5x2x2 ft.	6x2x2 ft.
External dimensions	2x3x7 ft.	3x3x7 ft.	4x3x7 ft.	5x3x7 ft.	6x3x7 ft.
Inflow velocity	105 fpm				
Down flow velocity	60 fpm				
Filter	HEPA filter 99.999% efficient				
Light	Fluorescent light				
Display	LED digital				
Body material	Wooden / SS 304 / SS 316				
Pressure gauge	Magnehelic gauge				

Blower Assembly	1/3 hp , 1440 rpm
Hepa filter	2no.
Castor wheels	Optional

www.oracleequipments.com