

LAMINAR AIR FLOW CABINET



Key specifications:-

Working area	2x2x2 ft.	3x2x2 ft.	4x2x2 ft.	6x2x2 ft.	8x2x2 ft.
HEPA filter size	2'x2'x6"	3'x2'x6"	4'x2'x6"	3'x2'x6"	4'x2'x6"
No. of hepa filters	1no.	1no.	1no.	2no.	2no
Air flow velocity	90 FPM				
HEPA filter	99.99% efficient				
Pre-filter	85% efficient				
Cabinet material	SS 304 / SS 316				
Light	UV c type				
Noise level	< 55 dba				
Power supply	220v , AC< single phase				

Castor wheels	Yes no.4
Blower	Yes no.1

www.oracleequipments.com