

HUMIDITY CHAMBER



Key specifications: -

Display	LED Digital
Power	220v , AC, Single Phase
Timer	
Capacity	95 Ltr
Accuracy	±1°C
Temperature Range	10-60°C
Humidity Range	50-90°C
Least Count	0.1°C
Body Material	MS Powder Coated

Highlights: -

No. Of specimen Keeping Racks	2no.
Wet Heater	1 Kw
Dry Heater	250 w
Temperature Sensor	PT 100 , 1no.
Humidity Sensor	1 (rh)
Outer Dimensions	
Inner Chamber Dimensions	18x18x18, 24x24x24, as per requirement
Finish	Blue Powder Coated
Castor wheels	4no.

Optional Features: -

- 1) Automatic Water filling.
- 2) Outer Body of Stainless Steel 304 grade.